



EZHold Lid Holder

The EZHold Lid Holders are accessories for the C.BIRD, C.NEST and S.NEST microbioreactors that keep lids from coming into direct contact with microwell plates. With the EZHold in place, you can easily add supplements or other components to wells while minimizing sample cross contamination.



Features

- Separates lid and microwell plate
- Keeps fluidic tubes from coming into contact
- Reduces the risk of cross contamination between wells
- Easy-to-use design
- Compatible with 96- and 24-well formats
- Available in C.BIRD, C.NEST and S.NEST

Specification

Dimension (WxDxH):
17.8 x 22.3 x 20 cm

Weight: 2.2 kg

EZHold ordering information

Compatibility	C.BIRD	C.NEST and S.NEST
Item No.	AC0104	AS0101

Learn more at: cytena-bps.com/ezhold/

CONTACT US

info@cytena-bps.com | (+886) 2-2720-6135 | www.cytene-bps.com

CYTENA BPS 
A BICO COMPANY